

The Open Stable Sling



Free your funny bone while providing comfort, stability and rotation control of the shoulder.

Top Shelf Orthopedics introduces an innovative approach to shoulder stabilization with a patent pending open elbow design.

The Open Stable Sling

Wearing a sling after a shoulder surgery used to be an uncomfortable experience, which meant less compliance and discomfort. By leaving the elbow open, these issues are addressed with **The Open Stable Sling**. Allowing the arm to sit in a natural but supported position helps reduce the chances of irritation to the nerves and soft tissue around the elbow while the shoulder heals.

Indications:

- Rotator Cuff Repair
- Bankart Repair
- Shoulder Instability
- Total Shoulder Replacement
- Capsule Repair
- Muscle / Ligament Injury
- SLAP Repair
- Shoulder Dislocation
- Posterior Capsule Repair
- AC Joint Separation

Features and Benefits:

- Universal lt/rt - reduces inventory
- Standard Abduction Pillow, External Rotation Abduction Pillow Options
- Internal, neutral & external rotation adjustable to accommodate rehab specific guidelines
- Multiple thumb support options to decrease irritation
- Comfortable shoulder strap reduces pain and pressure on the unaffected shoulder
- Trimmable strap ends allow for adjustment and optimal fit to the wearer
- Easy one handed buckle closure allows the wearer to don and doff the device without assistance
- Allows “Sling Only” configuration to aid in the recovery process when abduction is no longer needed



Ordering Information:

Description	Part Number	Size
The Open Stable Sling, Standard Abduction Pillow	502083	Small
The Open Stable Sling, Standard Abduction Pillow	502084	Medium
The Open Stable Sling, Standard Abduction Pillow	502085	Large
The Open Stable Sling, Standard Abduction Pillow	502086	X-Large
The Open Stable Sling, External Rotation Abduction Pillow	502043	Small
The Open Stable Sling, External Rotation Abduction Pillow	502044	Medium
The Open Stable Sling, External Rotation Abduction Pillow	502045	Large
The Open Stable Sling, External Rotation Abduction Pillow	502046	X-Large
The Open Stable Sling, Sling Only	502093	Small
The Open Stable Sling, Sling Only	502094	Medium
The Open Stable Sling, Sling Only	502095	Large
The Open Stable Sling, Sling Only	502096	X-Large
15 Degree Accessory Wedge	502060	N/A

Ph: 866.592.0488
Fax: 209.834.8832
1851 E. Paradise Rd., Ste A
Tracy, CA 95304

Our top priority is for patients to be safe, have positive outcomes, and recover as quickly as possible.

To view our full product line, please visit us at:

topshelfortho.com

© 2024 Top Shelf Orthopedics. All rights reserved. 05/2024 Rev. A

