


SIERRA COLD THERAPY SYSTEM™



Achieve the summit with the Sierra Cold Therapy System™

Introducing a new cold therapy system with customizable time and temperature settings. When it comes to cold therapy, not every patient is plug and play.

Our top priority is for patients to be safe, have positive outcomes, and recover as quickly as possible.

The Sierra Cold Therapy System™ follows these principles through proven technology and unique features that customize each therapy session.



The Sierra Cold Therapy System™ reduces edema, swelling, and pain where cold therapy is indicated for post-surgical and acute injuries.



Temperature Options:
45, 50 and 55 degrees



Timed Sessions:
20, 40 and 60 minutes



Quick Start Option:
One-touch to begin therapy



Additional Features:

- Automatic Shut off
- Lockable lid to reduce spillage
- Clear cooler with visible ice and water fill labels

Pad Options:

- Knee
- Shoulder
- Multi-use



Description	Part Number
Sierra Cold Therapy System Unit	508200
Sierra Cold Therapy System Unit w/ Multi-use Pad	508201
Sierra Cold Therapy System Unit w/ Shoulder Pad	508202
Shoulder Pad	508210
Multi-use Pad	508211
Knee Pad	508212
Power Supply	508220
Compression Straps for pads	508221
Power Supply, 3 Prong	508222



QR Code for access to instructions and troubleshooting.

The Sierra Cold Therapy System must be prescribed by a physician or licensed healthcare practitioner.

topshelfortho.com

© 2024 Top Shelf Orthopedics. All rights reserved. 06/2024 Rev. A

